



PRODUCT SAFETY DATA SHEET

Salts Healthcare Limited
Richard Street Aston Birmingham B7 4AA
Telephone: +44 (0) 121 333 2000
Fax: +44 (0) 121 333 2010

E-mail: salt@salts.co.uk Website: www.salts.co.uk

PRODUCT NAME(S):

Secuplast Hydro Aloe (SPHA2)

2. INTENDED USES:

As an additional adhesive for Ostomy pouches.

3. COMPOSITION:

A brown coloured putty like substance with Aloe Vera manufactured in strip form. Consisting of poly- iso butylene, gelatine, Aloe Vera and sodium carboxy-methyl-cellulose. Latex Free. Laminated to film silicone coated release liner. Top coating of PU film.

4. PHYSICAL AND CHEMICAL PROPERTIES:

Adhesive polymer, chemically inert.

5. HEALTH HAZARD:

No health hazard is anticipated during normal handling of this product. Avoid contact with eyes - may cause irritation. Stop using product if skin sensitisation occurs.

Tested for skin irritation (ISO 10993-10:1995):- Slight skin irritant

6. FIRE HAZARD AND EMERGENCY ACTION:

Flammable. Combustible solid which gives off toxic fumes, mainly carbon monoxide and carbon dioxide, when ignited. Extinguish fires with water, foam, dry chemicals, halons or carbon halons or carbon dioxide. Wear self-contained breathing apparatus if large quantities are ignited.







PRODUCT SAFETY DATA SHEET:-SECUPLAST HYDRO WITH ALOE

Page 2 of 2



7. STORAGE PRECAUTIONS:

Store in a cool dry place.

8. TRANSPORT PRECAUTIONS:

No special precautions are required.

9. HANDLING/USE/PROTECTIVE CLOTHING:

Minimal hazard.

10. DISPOSAL:

Discard used product safely. Controlled disposal of large quantities according to local environmental health guidelines.

11. FIRST AID:

a) Inhalation: Not applicableb) Contact with skin: Not applicable

c) Contact with eyes: Irrigate well with copious amounts of clean water. Seek medical aid.

d) Ingestion Wash out mouth with water. DO NOT induce vomiting. Seek medical aid.

12. ADDITIONAL INFORMATION:

Accredited by the British Skin Foundation

13. NAME, ADDRESS AND TELEPHONE NUMBER OF SUPPLIER

Salts Healthcare Ltd TEL: 0121 333 2000 Richard Street FAX: 0121 333 2020

Birmingham B7 4AA

14. REFERENCE NUMBER, DATE OF ISSUE:

SPHA2 October 2010

This information is provided in accordance with the requirements of the UK Health and Safety at Work Act 974 and specifically in order to assist risks' of the product to make their 'Assessment of health risks' as required by the UK Control of Substances Hazardous to Health Regulation 1988 (COSHH assessments). Provision of this information does not preclude users from seeking advice from other sources as indicated in the COSHH guides.

This information is intended to cover potential hazards at the place of work and does not detail medical uses, indications, contra-indications and precautions for the treatment of patients





