

PRODUCT SAFETY DATA SHEET





Salts Healthcare Limited
Richard Street Aston Birmingham B7 4AA
Telephone: +44 (0) 121 333 2000
Fax: +44 (0) 121 333 2010

E-mail: salt@salts.co.uk Website: www.salts.co.uk

PRODUCT NAME(S):

Elastic Belt

2. INTENDED USES:

For attachment to ostomy pouches

3. COMPOSITION:

An elastic fabric belt comprising 65% rayon,11% covered rubber Elastane and 24% cotton with polypropylene connectors and adjuster.

4. PHYSICAL AND CHEMICAL PROPERTIES:

Chemically inert.

5. HEALTH HAZARD:

No health hazard is anticipated during normal handling of this product.

6. FIRE HAZARD AND EMERGENCY ACTION:

Product will take fire.

7. STORAGE PRECAUTIONS:

No special storage requirements. Keep away from excessive heat sources.











PRODUCT SAFETY DATA SHEET – ELASTIC BELT Page 2 of 2 Uncontrolled Document



8. TRANSPORT PRECAUTIONS:

No special precautions.

9. HANDLING/USE/PROTECTIVE CLOTHING:

No special precautions.

10. DISPOSAL:

Controlled disposal according to local environmental health & waste disposal guidelines.

11. FIRST AID:

a) Inhalation:

Not applicable.

b) Contact with skin:

Not applicable.

c) Contact with eyes:

Not applicable.

d) Ingestion

Not applicable.

12. ADDITIONAL INFORMATION:

13. NAME, ADDRESS AND TELEPHONE NUMBER OF SUPPLIER

Salt & Son Ltd TEL: 0121 333 2000 Lord Street FAX: 0121 359 0830

Birmingham B7 4DS

14. REFERENCE NUMBER, DATE OF ISSUE:

AB01 JANUARY 2003

This information is provided in accordance with the requirements of the UK Health and Safety at Work Act 1974 and specifically in order to assist risks' of the product to make their 'Assessment of health risks' as required by the UK Control of Substances Hazardous to Health Regulation 1988 (COSHH assessments). Provision of this information does not preclude users from seeking advice from other sources as indicated in the COSHH guides.

This information is intended to cover potential hazards at the place of work and does not detail medical uses, indications, contra-indications and precautions for the treatment of patients. S:\Quality\COSHH\Belt AB01.doc

This Adobe Acrobat .pdf document is freely available on the Salts Server, it may be printed, but not altered. Please always check with QA that you are using the correct document.







