

PRODUCT DESCRIPTION

Nelaton Catheters for intermittent bladder drainage.

PRODUCT INFORMATION

- Sterile, individual peel pouches
- 100 pouches per inner polythene bag
- 38cm length
- Soft rounded tip for atraumatic insertion
- Two side eyes for maximum drainage
- Supplied single wrapped (SW) as standard
- Latex Free
- Classification: 1
- Material: Catheter = Medical grade PVC (Shore Hardness A80)

Single Wrapped Code No.	Qty Per Box	CH/Dia.
NC-1206/SW	500	6ch Ø 2.00mm
NC-1208/SW	500	8ch Ø 2.66mm
NC-1210/SW	500	10ch Ø 3.33mm
NC-1212/SW	500	12ch Ø 4.00mm
NC-1214/SW	500	14ch Ø 4.66mm
NC-1216/SW	500	16ch Ø 5.33mm
NC-1218/SW	500	18ch Ø 6.00mm
NC-1220/SW	400	20ch Ø 6.66mm
NC-1222/SW	400	22ch Ø 7.33mm
NC-1224/SW	400	24ch Ø 8.00mm

INTENDED USE AND DURATION OF USE

Intended use: Intermittent bladder drainage.

Duration of use: Intermittent but less than 60 minute duration.

SAFETY AND STORAGE INFORMATION

Toxicity: Manufactured from biocompatible materials.

Storage: Do not store below 3 degrees or above 30 degrees.

Important Notes:

The information given in this sheet is not intended to be exhaustive. The information contained in