
	TITLE: DEVICE MASTER RECORD	
	QUALITY FILE	Ref. no.: QF 04-5
	Revision: 3	Page 1 of 2
	Issue date: 15.07.1998	Effective date: 02.01.2023

SUCTION SET

REF. NO. 1113 2100 83 / 1113 3000 83 / 1113 3500 83



POOL SUCTION STRAIGHT

COMPONENT DESCRIPTION	COMPONENT NO.
1. SUCTION PIECE	
a. outer rigid tube ID/OD 8.0/10.0, length: 210 mm (20 holes, 2 holes and cut for easy connection)	0.08021010
b. inner rigid tube ID/OD 4.0/6.5, length: 335 mm, 4 relieving eyes	0.04033565
c. handle ID/OD 6.5/8.5, length: 125 mm, straight, PVC	0.06512585
3. TUBE	
Flexible transparent tube with ribs ID/OD 5.8/8.3 mm material: PVC 73°ShoreA, DEHP – FREE,	
• 1113 2100 83 – length: 1900 mm	0.58190083
• 1113 3000 83 – length: 2900 mm	0.58290083
• 1113 3500 83 – length: 3400 mm	0.58340083
4. PROXIMAL END	
funnel 8.3 material: PVC 50°ShoreA, DEHP – FREE, blue	0.00000083
PACKAGING MATERIAL	
Innerbag 165 x 550 mm	2.16555004
Pouch 210 x 490 mm	5.21049410
Label on the pouch 85 x 131 mm	7.00085131
Outerbag 310 x 600 mm	2.31060010
Box 386 x 386 x 260 mm	6.38638626
Label on the box 85 x 131 mm	7.00085131

25 pcs/bag

50 pcs/box



TITLE: DEVICE MASTER RECORD	
QUALITY FILE	Ref. no.: QF 04-5
Revision: 3	Page 2 of 2
Issue date: 15.07.1998	Effective date: 02.01.2023

Work instruction QI 01
Packaging instruction QI 04

Comments:	Glued	No glued
Outer rigid tube		X
Inner suction piece and tube	X	
Placing the product in an envelope as in the first photo.		