
World's strongest surgical blade? ... KIATO











Aditya Dispomed Products Pvt. Ltd.
Plot No. 19 Sector 6 IMT Manesar
Gurgaon Haryana 122050 India

Tel. +91 124 47 64 900 Fax +91 124 47 64 911

info@adityadispomed.com www.adityadispomed.com

A subsidiary of Kehr Surgical Pvt. Ltd.



Safe Back sterile surgical blades

Aditya Dispomed is pleased to offer a complete range of KIATO SAFE BACK Carbon Steel Surgical Blades featuring a patent pending safety innovative feature of a TWIN rib. These blades are on an average atleast 20% stronger than competing leading products. They can be used with a standard Bard Parker metal handle fitment 3 & 4 conforming to BS EN 227740.

A study of the complaints reported by users on MAUDE (FDA) site clearly shows that the incidence of blade breaking during use far exceeds all other incidences. We believe our innovative surgical blade suitably mitigates and reduces risk of breakage during use (surgery)

> Though statistically Aditya Dispomed Products Pvt. Ltd., reflects superior performance based on Maude data, we proactively (pre-priory) acted on above thread most prominent problems faced by uses about other devices available in international markets





Figure 4.1: Trend of adverse events reported on Madue (FDA) site regarding blade/ disposable scalpel year 2000-2016





Size

