sempermed®(s) derma PF

DERMA PF FOR SENSITIVE HANDS

HIGHLY SKIN-FRIENDLY – no powder-related allergies EASY TO DON – patented synthetic inner coating WORK CONFIDENTLY – micro-rough surface for great grip

SEMPERMED[®] DERMA PF

FOR SENSITIVE HANDS

The Sempermed[®] Derma PF offers good reasons to switch to powder-free gloves.

HIGHLY SKIN-FRIENDLY

This powder-free glove provides improved safety for staff, patients and visitors as there are no powder-related dermatological, allergic and toxic health risks.

EASY TO DON

The glove is easy to put on thanks to the synthetic lining.

SURE GRIP

The Sempermed[®] Derma PF offers a confident grip and thus combines safety and tactile properties.

AN ECONOMICALLY SOUND SOLUTION

Allergies and other health risks not only pose risks to end users, they can also put a significant strain on health care systems. This powder-free glove removes the risk of powderrelated health problems and thereby reduces potential costs in the long term.

TECHNICAL DATA

sempermed[®] () derma PF

ТҮРЕ	Sterile surgical glove for single use, powder-free, with synthetic lining
MATERIAL	Natural rubber latex (NR)
COLOUR	Natural white
GLOVE SHAPE	Fully anatomical with rolled rim
SIZE/OVERALL LENGTH as per EN 455-2	51/2, 6 and 61/2 min. 260 mm 7, 71/2 and 8 min. 270 mm 81/2 and 9 min. 280 mm
WALL THICKNESS measured in single layer in the palm area	0.16–0.20mm
BARRIER PERFORMANCE	AQL 0.65
FORCE AT BREAK exceeding EN 455-2 (≥ 12 N instead 9 N)	≥ 12 N
DURABILITY	3 years in original packaging if stored as per DIN 7716, ISO 2230
STERILISATION	Radiation STERILE R
LABELLING	EN 1041, EN 455, EN 420, EN 374, ISO 15223
PURPOSE according to MP Regulation 93/42/EWG according to PPE Regulation (EU) 2016/425	Medical device class IIa Single-use protective glove – PPE category III (protection against chemical substances for limited time) ¹⁾
PACKAGING	Left and right glove with turned up cuff in inner pouch, ozone-tight, sealed in flat medical fibre-free peel pack. In dispenser carton with sterilisation indicator: 50 pairs In transport carton with sterilisation indicator: 300 pairs
ARTICLE NUMBERS	Size 5 ½822751522Size 7 ½822751722Size 6822751602Size 8822751802Size 6 ½822751622Size 8½822751822Size 7822751702Size 9822751902

CONTACT US!

* Please note the Sempermed Chemical Resistance list is available on www.sempermed.com/sempermed_informs/chemical_resistance IMPORTANT NOTE: The latest product information is available at www.sempermed.com under the 'Products' tab. Failure to observe this information, in particular with regard to (chemical) resistance, frequency of use and tolerability of the gloves, can result in personal injury and/or material damage. Semperit is not liable for incorrect use of the gloves nor will Semperit be liable for improper storage and handling. In case of damage or unusual signs of wear and tear, gloves shall not be used. In case of any doubt, expert advice should be sought before use. The information and classification is up to date as of the date of publication. Technical details are average values from production and may vary in individual cases. This disclaimer is subject to mistakes, printing errors and amendments. CAUTION: Natural latex can cause allergic reactions, including anaphylactic shock. Semperit bears no liability for allergic reactions. [©]Copyright 2019 Sempermed