




Product data sheet - 60.01.01

24.11.2015

General Data			
Included in	<b>Catalogues</b> 450.00.00 (Kat. 12), 450.00.01 (HK), 450.00.02 (HK), 450.00.03 (HK neu), 450.00.25 (Kat. 1),		
Last update	2014-12-17		
Title	<b>Eye Speculum</b>		
Article number	<b>60.01.01</b>	Author	Barraquer
Description	Eye Speculum		
Length	4 cm - 1 1/2"	Weight	
		Diameter	
Properties	child size	Tip / Working end	
Article link		Sales unit	
Accessories		Number / Size	
Complete picture		Detail images	
			
Quality management data			
Preparation instructions	MPK01		
Operation manuals	<i>No operation manuals available.</i>		
GTIN	4046826096625		
GMDN	35349		
Risk category EC	1		
UMDNS	17-933		
Net weight	0.0020		



III CMF III SURGICAL III NEURO+SPINE


**medicon**<sup>®</sup>

Product data sheet - 60.01.01

24.11.2015

Materials			
	Material number	Description	Elements
Main component material	1.4310	X10CrNi18-8	C 0,05 - 0,15% Cr 16,0 - 19,0% Ni 6,0 - 9,5%
Turned part material			
Handle material			
Others parts material			
Coating material			

Delivery mode	
---------------	---