




Product data sheet - 76.01.03

24.11.2015

General Data			
Included in	<b>Catalogues</b> 450.00.01 (HK), 450.00.02 (HK), 450.00.03 (HK neu), 450.00.30 (Kat. 8), 450.00.47 (Cardio), 450.00.50 (Ortho), 450.00.52 (implantArt), 450.00.53 (Cardio neu), 450.00.54 (OMF Instr), 450.00.55 (OMF Impl), 450.00.57 (Plastic), 450.00.60 (Plastic neu), 450.00.61 (implantArt), 450.00.63 (Ortho), <b>Brochures</b> 451.09.01 (bimax), 451.09.02 (bimax),		
Last update	2014-08-07		
Title	<b>Bone Suturing Wire</b>		
Article number	<b>76.01.03</b>	Author	
Description	Bone Suturing Wire, <b>soft</b> , of austenitic stainless steel for implants, compatible with tissue		
Length		Weight	
		Diameter	Ø 0.30 mm - 29 G
Properties	coil/10 m - 33.3 ft.	Tip / Working end	
Article link		Sales unit	
Accessories		Number / Size	
Complete picture		Detail images	
			
Quality management data			
Preparation instructions	MPK01		
Operation manuals	G7080000_2014_04		
GTIN	4046826146542		
GMDN	35685		
Risk category EC	2b		
UMDNS	16-104		
Net weight	0.0170		



III CMF III SURGICAL III NEURO+SPINE


**medicon**<sup>®</sup>

Product data sheet - 76.01.03

24.11.2015

Materials			
	Material number	Description	Elements
Main component material	1.4441	X2CrNiMo18-15-3	Cr 17,0 - 18,5% Mo 2,5 - 3,2% Ni 13,0 - 15,5%
Turned part material			
Handle material			
Others parts material			
Coating material			

Delivery mode	
---------------	---