





Product data sheet - 60.05.13

24.11.2015

General Data			
Included in	Catalogues 450.00.00 (Kat. 12), 450.00.01 (HK), 450.00.02 (HK), 450.00.03 (HK neu), 450.00.30 (Kat. 8), 450.00.35 (Kat. 15), 450.00.37 (ENT neu), 450.00.47 (Cardio), 450.00.48 (Plastic), 450.00.53 (Cardio neu), 450.00.54 (OMF Instr), 450.00.57 (Plastic), 450.00.60 (Plastic neu),		
Last update	2008-01-25		
Title	Lid Retractor		
Article number	60.05.13	Author	Desmarres
Description	Lid Retractor		
Length	14 cm - 5 1/2"	Weight	
		Diameter	
Properties		Tip / Working end	10 x 16 mm
Article link		Sales unit	
Accessories		Number / Size	# 3
Complete picture		Detail images	
			
Quality management data			
Preparation instructions	MPK01		
Operation manuals	<i>No operation manuals available.</i>		
GTIN	4046826096991		
GMDN	46767		
Risk category EC	1		
UMDNS	13-381		
Net weight	0.0150		




III CMF III SURGICAL III NEURO+SPINE

medicon[®]

Product data sheet - 60.05.13

24.11.2015

Materials			
	Material number	Description	Elements
Main component material	1.4021	X20Cr13	C 0,16 - 0,25% Cr 12 - 14%
Turned part material			
Handle material			
Others parts material			
Coating material			

Delivery mode	
---------------	---