





Product data sheet - 60.74.03

24.11.2015

General Data			
Included in	Catalogues 450.00.00 (Kat. 12), 450.00.01 (HK), 450.00.02 (HK), 450.00.03 (HK neu), 450.00.30 (Kat. 8),		
Last update	2008-01-25		
Title	Micro Scissors		
Article number	60.74.03	Author	Mc Pherson Vannas
Description	Micro Scissors		
Length	9.5 cm - 3 3/4"	Weight	
		Diameter	
Properties	curved	Tip / Working end	sharp tips
Article link		Sales unit	
Accessories		Number / Size	
Complete picture		Detail images	
			
Quality management data			
Preparation instructions	MPK02		
Operation manuals	<i>No operation manuals available.</i>		
GTIN	4046826101008		
GMDN	13488		
Risk category EC	1		
UMDNS	13-488		
Net weight	0.0340		




III CMF III SURGICAL III NEURO+SPINE

medicon[®]

Product data sheet - 60.74.03

24.11.2015

Materials			
	Material number	Description	Elements
Main component material	1.4117	X38CrMoV15	C 0,35 - 0,40% Cr 14 - 15%
Turned part material	1.4305	X8CrNiS18-9	Cr 17,0 - 19,0% Ni 8,0 - 10,0%
Handle material			
Others parts material			
Coating material			

Delivery mode	
---------------	---