





Product data sheet - 04.33.70

24.11.2015

General Data			
Included in	Catalogues 450.00.00 (Kat. 12), 450.00.01 (HK), 450.00.02 (HK), 450.00.03 (HK neu), 450.00.35 (Kat. 15), 450.00.45 (Neuro), 450.00.47 (Cardio), 450.00.48 (Plastic), 450.00.53 (Cardio neu), 450.00.57 (Plastic), 450.00.58 (Cranialis), 450.00.59 (Spinalis), 450.00.60 (Plastic neu), Brochures 451.03.00 (Plastisch), 451.03.23 (Medicut),		
Last update	2008-01-24		
Title	Dissecting Scissors		
Article number	04.33.70	Author	Metzenbaum "HM"
Description	Dissecting Scissors		
Length	20 cm - 8"	Weight	
		Diameter	
Properties	curved	Tip / Working end	
Article link		Sales unit	
Accessories		Number / Size	
Complete picture		Detail images	
			
Quality management data			
Preparation instructions	MPK02		
Operation manuals	<i>No operation manuals available.</i>		
GTIN	4046826004200		
GMDN	13492		
Risk category EC	1		
UMDNS	13-492		
Net weight	0.0560		



III CMF III SURGICAL III NEURO+SPINE

medicon[®]

Product data sheet - 04.33.70

24.11.2015

Materials			
	Material number	Description	Elements
Main component material	1.4117	X38CrMoV15	C 0,35 - 0,40% Cr 14 - 15%
Turned part material	1.4305	X8CrNiS18-9	Cr 17,0 - 19,0% Ni 8,0 - 10,0%
Handle material			
Others parts material	WSG_20-GO-55_CRTZ	RCoCr-C	
Coating material	T_GOLD	Teil des Instrumentes ist Vergoldet Goldschicht galvanisch max 0,004mm dick	

Delivery mode

