





Product data sheet - 204.33.14

24.11.2015

General Data			
Included in	Catalogues 450.00.01 (HK), 450.00.02 (HK), 450.00.03 (HK neu), 450.00.30 (Kat. 8), 450.00.35 (Kat. 15), 450.00.36 (ENT), 450.00.37 (ENT neu), 450.00.45 (Neuro), 450.00.47 (Cardio), 450.00.48 (Plastic), 450.00.53 (Cardio neu), 450.00.57 (Plastic), 450.00.58 (Cranialis), 450.00.59 (Spinalis), 450.00.60 (Plastic neu), Brochures 451.03.23 (Medicut),		
Last update	2014-10-16		
Title	Dissecting Scissors		
Article number	204.33.14	Author	Metzenbaum supercut
Description	Operating and Dissecting Scissors		
Length	14 cm - 5 1/2"	Weight	
		Diameter	
Properties	curved	Tip / Working end	
Article link		Sales unit	
Accessories		Number / Size	
Complete picture		Detail images	
			
Quality management data			
Preparation instructions	MPK02		
Operation manuals	<i>No operation manuals available.</i>		
GTIN	4046826172398		
GMDN	13492		
Risk category EC	1		
UMDNS	13-492		
Net weight	0.0300		




III CMF III SURGICAL III NEURO+SPINE

medicon[®]

Product data sheet - 204.33.14

24.11.2015

Materials			
	Material number	Description	Elements
Main component material	1.4117	X38CrMoV15	C 0,35 - 0,40% Cr 14 - 15%
Turned part material	1.4305	X8CrNiS18-9	Cr 17,0 - 19,0% Ni 8,0 - 10,0%
Handle material			
Others parts material			
Coating material	T_GOLD	Teil des Instrumentes ist Vergoldet Goldschicht galvanisch max 0,004mm dick	

Delivery mode	
---------------	--