




Product data sheet - 45.51.22

24.11.2015

General Data			
Included in	<b>Catalogues</b> 450.00.00 (Kat. 12), 450.00.01 (HK), 450.00.02 (HK), 450.00.03 (HK neu), 450.00.35 (Kat. 15), 450.00.36 (ENT), 450.00.37 (ENT neu),		
Last update	2011-02-01		
Title	<b>Frame, only</b>		
Article number	<b>45.51.22</b>	Author	Davis-Boyle
Description	Frame only		
Length		Weight	
		Diameter	
Properties	standard size, All available frames and tongue blades for the mouth gags can be combined according to choice.	Tip / Working end	
Article link		Sales unit	
Accessories		Number / Size	
Complete picture		Detail images	
			
Quality management data			
Preparation instructions	MPK01		
Operation manuals	<i>No operation manuals available.</i>		
GTIN	4046826048976		
GMDN	35085		
Risk category EC	1		
UMDNS	11-822		
Net weight	0.1070		



III CMF III SURGICAL III NEURO+SPINE

**medicon**<sup>®</sup>

Product data sheet - 45.51.22

24.11.2015

Materials			
	Material number	Description	Elements
Main component material	1.4021	X20Cr13	C 0,16 - 0,25% Cr 12 - 14%
Turned part material			
Handle material	1.4301	X5CrNi18-10	Cr 18,0 - 19,5% Ni 8,0 - 10,5%
Others parts material	1.4310	X10CrNi18-8	C 0,05 - 0,15% Cr 16,0 - 19,0% Ni 6,0 - 9,5%
Coating material			

Delivery mode

