





Product data sheet - 09.08.13

24.11.2015

General Data			
Included in	Catalogues 450.00.00 (Kat. 12), 450.00.01 (HK), 450.00.02 (HK), 450.00.03 (HK neu), 450.00.35 (Kat. 15), 450.00.36 (ENT), 450.00.37 (ENT neu), 450.00.40 (Kat. 5), 450.00.54 (OMF Instr),		
Last update	2008-01-23		
Title	Myrtle leaf Probe		
Article number	09.08.13	Author	
Description	Myrtle leaf Probe		
Length	13 cm - 5 1/8	Weight	
		Diameter	
Properties		Tip / Working end	
Article link		Sales unit	
Accessories		Number / Size	
Complete picture		Detail images	
			
Quality management data			
Preparation instructions	MPK01		
Operation manuals	<i>No operation manuals available.</i>		
GTIN	4046826010478		
GMDN	33433		
Risk category EC	1		
UMDNS	13-117		
Net weight	0.0040		




III CMF III SURGICAL III NEURO+SPINE

medicon[®]

Product data sheet - 09.08.13

24.11.2015

Materials			
	Material number	Description	Elements
Main component material	1.4305	X8CrNiS18-9	Cr 17,0 - 19,0% Ni 8,0 - 10,0%
Turned part material			
Handle material			
Others parts material			
Coating material			

Delivery mode	
---------------	---