





Product data sheet - 66.03.51

24.11.2015

General Data			
Included in	Catalogues 450.00.00 (Kat. 12), 450.00.35 (Kat. 15), 450.00.36 (ENT), 450.00.37 (ENT neu),		
Last update	2012-03-19		
Title	Nasal Speculum		
Article number	66.03.51	Author	Hartmann-Halle
Description	Nasal Speculum		
Length	15 cm - 6	Weight	
		Diameter	
Properties	aseptic joint (dismantable)	Tip / Working end	
Article link		Sales unit	
Accessories		Number / Size	# 1
Complete picture		Detail images	
			
Quality management data			
Preparation instructions	MPK01		
Operation manuals	<i>No operation manuals available.</i>		
GTIN	4046826109875		
GMDN	35350		
Risk category EC	1		
UMDNS	13-664		
Net weight	0.0770		




III CMF III SURGICAL III NEURO+SPINE

medicon[®]

Product data sheet - 66.03.51

24.11.2015

Materials			
	Material number	Description	Elements
Main component material	1.4021	X20Cr13	C 0,16 - 0,25% Cr 12 - 14%
Turned part material			
Handle material			
Others parts material			
Coating material			

Delivery mode	
---------------	---