





Product data sheet - 10.60.64

24.11.2015

General Data			
Included in	Catalogues 450.00.00 (Kat. 12), 450.00.01 (HK), 450.00.02 (HK), 450.00.03 (HK neu), 450.00.30 (Kat. 8), 450.00.35 (Kat. 15), 450.00.36 (ENT), 450.00.37 (ENT neu), 450.00.40 (Kat. 5),		
Last update	2008-01-25		
Title	Needle Holder		
Article number	10.60.64	Author	Olsen-Hegar "HM"
Description	Needle Holder and Scissors, standard pattern, with Tungsten Carbide inserts, gold-plated ring handles		
Length	14 cm - 5 1/2	Weight	
		Diameter	
Properties	fine cross serrated, straight	Tip / Working end	
Article link		Sales unit	
Accessories		Number / Size	
Complete picture		Detail images	
			
Quality management data			
Preparation instructions	MPK01		
Operation manuals	<i>No operation manuals available.</i>		
GTIN	4046826017163		
GMDN	12726		
Risk category EC	1		
UMDNS	12-726		
Net weight	0.0360		



III CMF III SURGICAL III NEURO+SPINE

medicon[®]

Product data sheet - 10.60.64

24.11.2015

Materials			
	Material number	Description	Elements
Main component material	1.4021	X20Cr13	C 0,16 - 0,25% Cr 12 - 14%
Turned part material			
Handle material			
Others parts material	K40	WcCo	
Coating material	T_GOLD	Teil des Instrumentes ist Vergoldet Goldschicht galvanisch max 0,004mm dick	

Delivery mode

