





Product data sheet - 10.60.67

24.11.2015

General Data			
Included in	Catalogues 450.00.00 (Kat. 12), 450.00.01 (HK), 450.00.02 (HK), 450.00.03 (HK neu), 450.00.35 (Kat. 15), 450.00.36 (ENT), 450.00.37 (ENT neu), 450.00.40 (Kat. 5),		
Last update	2008-01-25		
Title	Needle Holder		
Article number	10.60.67	Author	Olsen-Hegar "HM"
Description	Combined Needle Holder and Scissors, standard pattern, with Tungsten Carbide inserts, gold-plated ring handles		
Length	17 cm - 6 3/4	Weight	
		Diameter	
Properties	fein cross serrated, straight	Tip / Working end	
Article link		Sales unit	
Accessories		Number / Size	
Complete picture		Detail images	
			
Quality management data			
Preparation instructions	MPK02		
Operation manuals	<i>No operation manuals available.</i>		
GTIN	4046826017170		
GMDN	12726		
Risk category EC	1		
UMDNS	12-726		
Net weight	0.0460		




III CMF III SURGICAL III NEURO+SPINE

medicon[®]

Product data sheet - 10.60.67

24.11.2015

Materials			
	Material number	Description	Elements
Main component material	1.4021	X20Cr13	C 0,16 - 0,25% Cr 12 - 14%
Turned part material			
Handle material			
Others parts material	K40	WcCo	
Coating material	T_GOLD	Teil des Instrumentes ist Vergoldet Goldschicht galvanisch max 0,004mm dick	

Delivery mode	
---------------	--