





Product data sheet - 06.02.11

24.11.2015

General Data			
Included in	Catalogues 450.00.00 (Kat. 12), 450.00.01 (HK), 450.00.02 (HK), 450.00.03 (HK neu), 450.00.25 (Kat. 1), 450.00.35 (Kat. 15), 450.00.36 (ENT), 450.00.37 (ENT neu), 450.00.46 (Sinus), 450.00.54 (OMF Instr),		
Last update	2008-01-24		
Title	Dressing Forceps		
Article number	06.02.11	Author	
Description	Dressing Forceps		
Length	11.5 cm - 4 1/2	Weight	
		Diameter	
Properties	straight, slim	Tip / Working end	
Article link		Sales unit	
Accessories		Number / Size	
Complete picture		Detail images	
			
Quality management data			
Preparation instructions	MPK01		
Operation manuals	<i>No operation manuals available.</i>		
GTIN	4046826006235		
GMDN	11777		
Risk category EC	1		
UMDNS	14-257		
Net weight	0.0140		




III CMF III SURGICAL III NEURO+SPINE

medicon[®]

Product data sheet - 06.02.11

24.11.2015

Materials			
	Material number	Description	Elements
Main component material	1.4024	X15Cr13	C 0,12 - 0,17% Cr 12 - 14%
Turned part material			
Handle material			
Others parts material			
Coating material			

Delivery mode	
---------------	---