





Product data sheet - 06.05.11

24.11.2015

General Data			
Included in	Catalogues 450.00.00 (Kat. 12), 450.00.01 (HK), 450.00.02 (HK), 450.00.03 (HK neu), 450.00.25 (Kat. 1), 450.00.35 (Kat. 15), 450.00.36 (ENT), 450.00.37 (ENT neu),		
Last update	2008-01-24		
Title	Tissue Forceps		
Article number	06.05.11	Author	Standard
Description	Tissue Forceps		
Length	11.5 cm - 4 1/2	Weight	
		Diameter	
Properties	straight	Tip / Working end	1 x 2 teeth
Article link		Sales unit	
Accessories		Number / Size	
Complete picture		Detail images	
			
Quality management data			
Preparation instructions	MPK01		
Operation manuals	<i>No operation manuals available.</i>		
GTIN	4046826006433		
GMDN	11785		
Risk category EC	1		
UMDNS	14-257		
Net weight	0.0170		




III CMF III SURGICAL III NEURO+SPINE

medicon[®]

Product data sheet - 06.05.11

24.11.2015

Materials			
	Material number	Description	Elements
Main component material	1.4024	X15Cr13	C 0,12 - 0,17% Cr 12 - 14%
Turned part material			
Handle material			
Others parts material			
Coating material			

Delivery mode	
---------------	---