





Product data sheet - 240.25.20

24.11.2015

General Data			
Included in	Catalogues 450.00.01 (HK), 450.00.02 (HK), 450.00.03 (HK neu), 450.00.30 (Kat. 8), 450.00.50 (Ortho), 450.00.63 (Ortho), Brochures 451.03.23 (Medicut),		
Last update	2014-04-28		
Title	Bandage Scissors		
Article number	240.25.20	Author	Lister supercut
Description	Bandage Scissors		
Length	20 cm - 8	Weight	
		Diameter	
Properties		Tip / Working end	angled to side
Article link		Sales unit	
Accessories		Number / Size	
Complete picture		Detail images	
			
Quality management data			
Preparation instructions	MPK02		
Operation manuals	<i>No operation manuals available.</i>		
GTIN	4046826175221		
GMDN	13481		
Risk category EC	1		
UMDNS	13-481		
Net weight	0.1210		



III CMF III SURGICAL III NEURO+SPINE

medicon[®]

Product data sheet - 240.25.20

24.11.2015

Materials			
	Material number	Description	Elements
Main component material	1.4117	X38CrMoV15	C 0,35 - 0,40% Cr 14 - 15%
Turned part material	1.4305	X8CrNiS18-9	Cr 17,0 - 19,0% Ni 8,0 - 10,0%
Handle material			
Others parts material			
Coating material	T_GOLD	Teil des Instrumentes ist Vergoldet Goldschicht galvanisch max 0,004mm dick	

Delivery mode

