





Product data sheet - 29.51.52

24.11.2015

General Data			
Included in	<b>Catalogues</b> 450.00.00 (Kat. 12), 450.00.01 (HK), 450.00.02 (HK), 450.00.03 (HK neu), 450.00.25 (Kat. 1), 450.00.30 (Kat. 8), 450.00.35 (Kat. 15), 450.00.36 (ENT), 450.00.37 (ENT neu), 450.00.40 (Kat. 5), 450.00.50 (Ortho), 450.00.54 (OMF Instr), 450.00.63 (Ortho),		
Last update	2012-04-17		
Title	<b>Sharp Curette, oval</b>		
Article number	<b>29.51.52</b>	Author	Volkmann
Description	Sharp Curette, oval		
Length	17 cm - 6 3/4	Weight	
		Diameter	
Properties		Tip / Working end	
Article link		Sales unit	
Accessories		Number / Size	# 2
Complete picture		Detail images	
			
Quality management data			
Preparation instructions	MPK01		
Operation manuals	<i>No operation manuals available.</i>		
GTIN	4046826041380		
GMDN	31335		
Risk category EC	1		
UMDNS	11-084		
Net weight	0.0530		



III CMF III SURGICAL III NEURO+SPINE

**medicon**<sup>®</sup>

Product data sheet - 29.51.52

24.11.2015

Materials			
	Material number	Description	Elements
Main component material	1.4021	X20Cr13	C 0,16 - 0,25% Cr 12 - 14%
Turned part material			
Handle material	1.4301	X5CrNi18-10	Cr 18,0 - 19,5% Ni 8,0 - 10,5%
Others parts material			
Coating material			

Delivery mode

