





III CMF III SURGICAL III NEURO+SPINE

medicon[®]

Product data sheet - 62.14.01

24.11.2015

General Data			
Included in	Catalogues 450.00.00 (Kat. 12), 450.00.35 (Kat. 15), 450.00.36 (ENT), 450.00.37 (ENT neu),		
Last update	2011-03-17		
Title	Ear loop		
Article number	62.14.01	Author	Billeau
Description	Ear loop		
Length	15,5 cm - 6 1/8"	Weight	
		Diameter	
Properties		Tip / Working end	
Article link		Sales unit	
Accessories		Number / Size	# 1
Complete picture		Detail images	
			
Quality management data			
Preparation instructions	MPK01		
Operation manuals	<i>No operation manuals available.</i>		
GTIN	4046826102685		
GMDN	31020		
Risk category EC	1		
UMDNS	13-631		
Net weight	0.0110		



III CMF III SURGICAL III NEURO+SPINE

medicon[®]

Product data sheet - 62.14.01

24.11.2015

Materials			
	Material number	Description	Elements
Main component material	1.4401	X5CrNiMo17-12-2	Cr 16,5 - 18,5% Ni 10,0 - 13,0% Mo 2,0 - 2,5%
Turned part material			
Handle material	1.4104	X14CrMoS17	C 0,10 - 0,17% Cr 15,5 - 17,5% Mo 0,20 - 0,60%
Others parts material			
Coating material			

Delivery mode

