





Product data sheet - 62.35.12

24.11.2015

General Data			
Included in	<b>Catalogues</b> 450.00.00 (Kat. 12), 450.00.35 (Kat. 15), 450.00.36 (ENT), 450.00.37 (ENT neu),		
Last update	2014-01-30		
Title	<b>MICRO Ear Forceps</b>		
Article number	<b>62.35.12</b>	Author	
Description	MICRO Ear Forceps		
Length	8 cm - 3 1/8"	Weight	
		Diameter	
Properties	extremely delicate, dismantable shaft	Tip / Working end	4 x 0.4 mm, serrated jaws
Article link		Sales unit	
Accessories		Number / Size	
Complete picture		Detail images	
			
Quality management data			
Preparation instructions	MPK02		
Operation manuals	<i>No operation manuals available.</i>		
GTIN	4046826104665		
GMDN	35802		
Risk category EC	1		
UMDNS	16-451		
Net weight	0.0160		



III CMF III SURGICAL III NEURO+SPINE

**medicon**<sup>®</sup>

Product data sheet - 62.35.12

24.11.2015

Materials			
	Material number	Description	Elements
Main component material	1.4021	X20Cr13	C 0,16 - 0,25% Cr 12 - 14%
Turned part material			
Handle material	1.4021	X20Cr13	C 0,16 - 0,25% Cr 12 - 14%
Others parts material	1.4310	X10CrNi18-8	C 0,05 - 0,15% Cr 16,0 - 19,0% Ni 6,0 - 9,5%
Coating material			

Delivery mode

