




Product data sheet - 01.27.03

24.11.2015

General Data			
Included in	Catalogues 450.00.00 (Kat. 12), 450.00.01 (HK), 450.00.02 (HK), 450.00.03 (HK neu), 450.00.25 (Kat. 1), 450.00.30 (Kat. 8), 450.00.35 (Kat. 15), 450.00.36 (ENT), 450.00.37 (ENT neu), 450.00.40 (Kat. 5), 450.00.45 (Neuro), 450.00.48 (Plastic), 450.00.54 (OMF Instr), 450.00.57 (Plastic), 450.00.58 (Cranialis), 450.00.59 (Spinalis), 450.00.60 (Plastic neu), Brochures 451.03.00 (Plastisch), 451.09.01 (bimax), 451.09.02 (bimax), Single sheets 451.06.90 (Cristaca),		
Last update	2013-08-16		
Title	Scalpel handle		
Article number	01.27.03	Author	
Description	Scalpel Handle for exchangeable scalpel blades fig. 10, 11, 12, 12d, 15, 15c		
Length	12.5 cm - 5"	Weight	
		Diameter	
Properties		Tip / Working end	
Article link		Sales unit	
Accessories		Number / Size	no. 3
Complete picture		Detail images	
			
Quality management data			
Preparation instructions	MPK01		
Operation manuals	<i>No operation manuals available.</i>		
GTIN	4046826000165		
GMDN	12235		
Risk category EC	1		
UMDNS	12-236		
Net weight	0.0210		




III CMF III SURGICAL III NEURO+SPINE

medicon[®]

Product data sheet - 01.27.03

24.11.2015

Materials			
	Material number	Description	Elements
Main component material	1.4016	X6Cr17	C
Turned part material			
Handle material			
Others parts material			
Coating material			

Delivery mode	
---------------	---