





Product data sheet - 03.06.16

24.11.2015

General Data			
Included in	Catalogues 450.00.00 (Kat. 12), 450.00.01 (HK), 450.00.02 (HK), 450.00.03 (HK neu), 450.00.25 (Kat. 1), 450.00.35 (Kat. 15), 450.00.40 (Kat. 5), 450.00.45 (Neuro), 450.00.58 (Cranialis), 450.00.59 (Spinalis),		
Last update	2008-01-24		
Title	Operating Scissors		
Article number	03.06.16	Author	standard
Description	Operating Scissors		
Length	16.5 cm - 6 1/2	Weight	
		Diameter	
Properties	straight	Tip / Working end	sharp/blunt
Article link		Sales unit	
Accessories		Number / Size	
Complete picture		Detail images	
			
Quality management data			
Preparation instructions	MPK02		
Operation manuals	<i>No operation manuals available.</i>		
GTIN	4046826002671		
GMDN	38727		
Risk category EC	1		
UMDNS	13-480		
Net weight	0.0610		



III CMF III SURGICAL III NEURO+SPINE

medicon[®]

Product data sheet - 03.06.16

24.11.2015

Materials			
	Material number	Description	Elements
Main component material	1.4117	X38CrMoV15	C 0,35 - 0,40% Cr 14 - 15%
Turned part material	1.4305	X8CrNiS18-9	Cr 17,0 - 19,0% Ni 8,0 - 10,0%
Handle material			
Others parts material			
Coating material			

Delivery mode

