

IMMEDIATE HEMOGLOBIN RESULTS FROM A HAND-HELD ANALYZER



Results in less than 2 seconds

Reagent-free microcuvettes

40 days or 10,000 tests continuous use

No start-up process



FAST AND ACCURATE

- Results in less than 2 seconds
- Measuring range: 0.3 – 25.5 g/dL
- Imprecision: CV < 1%
- Factory calibrated against HiCN reference method according to ICSH
- Photometric method with broad spectrum and compensation for turbidity and scattering
- No re-calibration required

EASY TO USE

- 'Always on' technology. The analyzer is always ready for use
- Automatic self-check between every measurement
- Long microcuvette shelf life: 2.5 years from date of production
- Microcuvettes unaffected by moisture or temperature and can be used until expiration date after opening the bag

PRACTICAL AND PORTABLE

- Sample volume: <10 µL
- Capillary or venous blood (EDTA and Heparin anticoagulant can be used)
- Minimal maintenance
- Operating ambient temperature: 10° – 35° C
- Storage temperature: 0° – 50° C (short term -30° C – 70° C)
- Integrated rechargeable lithium-ion battery (40 days / 10,000 tests of continuous use)
- Compact: 150 (d) x 90 (w) x 40mm (h)
- Lightweight: 185 g
- Connect to PC via USB 2.0 or Bluetooth (optional)

Fast, easy to use and reliable in all climate conditions



1

Collect blood sample



2

Ensure analyzer is ready for use



3

Results appear in 1 - 2 seconds

Distributed by

Manufacturer

DiaSpect Medical GmbH
 Von-Cancrin-Str. 3
 63877 Sailauf
 Germany
Tel. +49 6093 996689-60

Sales

EKF Diagnostics Holdings plc
 Avon House, 19 Stanwell Road,
 Penarth, Cardiff, UK, CF64 2EZ
Tel. +44 (0)29 20 710570
www.ekfdiagnostics.com

Approved for the European Market

