## **EQUIPMENT DRAPES**

A range of sterile drapes for scopes and other imaging systems



## HIGH-DEFINITION CAMERA DRAPE

### The next level of camera drapes.

- Suitable for high-definition camera systems without compromising the HD image quality.
- Closed-end offers camera full protection from fluids and acts as a reliable sterile barrier.
- Allows quick and easy scope changes during procedures without compromising sterility.
- Broad-edged ring ensures ease of use and prevents contamination of the sterile field.
- · Sterile, single-use.

CODE	DIMENSIONS	UNIT OF SALE
ECN150-250CAP	Size 150mm x 2.5m, Ring aperture 85mm	50
TRAINING VIDEOS AVAILABLE		
*		





## 01 Attach the endoscope

Simply push the endoscope into the drape holder.





## 02 Attach the camera

Connect the camera onto the camera drape holder and extend the drape over the cable by pulling back the ring.

Any excess drape is retained inside the ring.

## QUICK CHANGE CAMERA DRAPE

- · Simple and secure 2-step fixing method.
- Allows quick and easy scope changes during procedures without compromising sterility.
- Closed-end offers cameras full protection from fluids and acts as a reliable sterile barrier.
- Broad-edged ring reduces risk of contaminating the sterile field.
- · Sterile, single-use.

CODE	DIMENSIONS	UNIT OF SALE
ECE150-250C	Size 150mm x 2.5m, Ring aperture 85mm	50





## 01 Attach the camera

Simple push-fit. The drape can then be extended over camera and cable by simply pulling back the ring.

Any excess drape is retained inside the ring.



02 Attach the endoscope

Simple and secure snap—on fitting.



## CONCENTRIC RING STYLE

- Ring style drapes ensure easy and efficient application without compromising sterility.
- Drape is fully contained within the ring and can be simply pulled back over camera and cable.
- · Any excess material remains securely inside the ring.
- Available with 3 different closures.
- · Sterile, single-use.

TRAINING VIDEOS AVAILABLE	

CODE	DIMENSIONS	UNIT OF SALE
ECN150-250	With two adhesive tape closures Size 150mm x 2.5m, Ring aperture 85mm	50
ECN150-250S	Elastic collar with two twist-lock spring clips and one adhesive tape closure Size 150mm x 2.5m, Ring aperture 85mm	50
ECN150-250TK	Elastomeric end section with two adhesive tape closures Size 150mm x 2.5m, Ring aperture 85mm	50



## **CONCERTINA STYLE**

- · Locator insert assists hand access through the drape.
- Folded cuff to facilitate application without compromising sterility
- · Sterile, single-use.

CODE	DIMENSIONS	UNIT OF SALE
EOD150-300	Supplied with two adhesive tape closures Size 150mm x 3m	25
EOD165-300	Supplied with two adhesive tape closures Size 165mm x 3m	25



## TELESCOPIC STYLE

- · Folded cuff to facilitate application without compromising sterility.
- · Sterile, single-use.

CODE	DIMENSIONS	UNIT OF SALE
EOD130-250T	With perforated tip and two adhesive tape closures Size 130mm x 2.5m	25
EOD165-250TK	With elastomeric end section and one adhesive tape closure Size 165mm x 2.5m, Aperture 10mm	25
EOD130-250TB	With cardboard sleeve Size 130mm x 2.5m	25







## FLAT FOLDED STYLE

- Folded cuff to facilitate application without compromising sterility
- · Sterile, single-use.

CODE	DIMENSIONS	UNIT OF SALE
EOD130-250	Supplied with one adhesive tape closure Size 130mm x 2.5m	25



## DRILL DRAPE - TELESCOPIC STYLE

- Folded cuff to facilitate application without compromising sterility.
- · Sterile, single-use.

CODE	DIMENSIONS	UNIT OF SALE
EOD130-100T	Supplied with two adhesive tape closure Size 130mm x 1m	20



## LAPAROSCOPIC INSTRUMENT POUCH

- Specially designed sterile pouch with five pockets to provide temporary storage for instruments during surgical procedures.
- Sterile, single-use.

CODE	DIMENSIONS	UNIT OF SALE
LIP500-048	With adhesive area to attach pouch to drapes/surfaces Size 500 x 480mm	25





# P<sup>3</sup>

## **OBJECTIVE LENS COVER**

- With removable lens protector.
- For use in conjunction with the P3 universal fitting microscope drape.
- · Sterile, single-use.

CODE	DIMENSIONS	UNIT OF SALE
MLC 048	48mm dia (for Zeiss OPMI)	10
MLC 065	65mm dia (for Zeiss MD)	10
MLC 070	70mm dia (for LEICA /WILD)	10



## MICROSCOPE HANDLE DRAPES

- · Clear plastic with adhesive fixing tape.
- · Packed in pairs.
- · Sterile, single-use.

CODE	DIMENSIONS	UNIT OF SALE
MHD180-330	Size 180 x 330mm	50 Pairs



## ULTRASOUND PROBE DRAPES

- Reduces the risk of patient-to-patient transmission of infectious agents.
- · Sterile, single-use.



## TELESCOPIC STYLE

CODE	DIMENSIONS	UNIT OF SALE
ULT127-150T	With two adhesive tape closures Size 127mm x 1.5m	25
ULT127-150T-G	Supplied with sterile gel sachets Size 127mm x 1.5m	25



## FLAT FOLDED STYLE

CODE	DIMENSIONS	UNIT OF SALE
UGD130-082	With two adhesive tape closuress Size 130 x 820mm	25



# P<sup>3</sup>

## **BONNET/BANDED BAG**

- Elasticated collar for easy application and increased security.
- With many uses in operating theatres and X-ray.
- · Sterile, single-use.

CODE	DIMENSIONS UNIT	IT OF SALE
CAB070	700mm dia 25	
CAB080	800mm dia 25	
CAB090	900mm dia 25	
CAB100	1.0m dia 25	
	P3 (E020)  See the second of t	


# OTHER EQUIPMENT DRAPES



## X-RAY SCREEN COVER

- Transparent film to allow clear view of screen during sterile procedures.
- Supplied with four adhesive tapes (two full-width) for a secure fit.
- · Sterile, single-use.


CODE	DIMENSIONS	UNIT OF SALE
SCC090-120	Size 900 x 1200mm	25



## MEMBRANE KEYPAD DRAPE

- Simple clear cover for enhanced security when using equipment keypad during sterile procedures.
- · Supplied with four adhesive fixing tapes.
- · Sterile, single-use.

CODE	DIMENSIONS	UNIT OF SALE
KPD200-210	Size 200 x 210mm	25



## COLPOSCOPE DRAPE, UNIVERSAL FITTING

- Clear plastic with elastomeric panel (140mm o/dia, 40mm i/dia) to fit more than one objective lens size.
- Single-use, packed under cleanroom conditions

CODE	DIMENSIONS	UNIT OF SALE
ID0001NS	Size 840 x 760mm	100

