



To whom it may concern

Stability of Visirub


We herewith confirm that Visirub has a shelf life of 1 year.


The ready-to-use solution is stable for three months after opening when stored in a dry and dark place. Partially used test preparation can be used again within this period of time.

Our statement is only valid for the original product container. Filling or refilling from or into other containers may not be carried out due to risk of contamination and intermixture.

Hamburg, February 12, 2010

BODE CHEMIE GmbH


Dr. Gerhard Bredthauer
Quality Control


Dr. Sven Eggerstedt
Development