



SURECLIP™


WITH OPENING SPANS OF 16, 11 AND 8 MM

The innovative SureClips™ from MICRO-TECH combine precise positioning with a reliable holding force. Therefore, the system is equally well suited for effective hemostasis as it is for endoscopic marking. Our clips can be rotated 360° and repeatedly open and close, which facilitates an enormously exact repositioning of the clips. Three versions are available: The SureClip™ Plus model with a span of 16 mm and fenestrated jaws and two other versions with 11 or 8 mm wingspan. Another important advantage of the

SureClip™ is the length of the remaining segment – the so called stem once the system has been withdrawn. Compared to all other available clip systems, the remaining segment with only 5 or 6 mm is up to 60% shorter than other clips. This is particularly advantageous in cases where multiple clips are being applied since it allows for a clearer endoscopic view during the operation and minimizes patient discomfort.

SPECIFIC CHARACTERISTICS

- 360° rotation for precise placement
- Maximum span of 16 mm
- Exclusively unlimited opening and closing cycles
- Powerful holding force
- Smooth handling, also in cases of difficult anatomical areas
- Ease of operation



Short stem lengths




The clips allow reopening and closing



COLORED HANDLES FOR EASY RECOGNITION

To recognize the different sizes SureClip™ quickly and easily, each size has uniquely colored handles and coating.



ECO-CLIP AS A COST-SAVING ALTERNATIVE

The uncovered Eco-Clips with a one-ring handle, left-handonly rotation, reclosable and repositionable are a cost-saving alternative to the SureClips™ and multi-clip systems.



PREMIUM ECO CLIP

The Premium Eco Clip offers the best results at the most economical conditions. The clips have a span of 11 or 16 mm. The short stem length helps placement in narrow lumens.



SPECIFICATIONS

REF	Clip opened mm	Length	Clip angle	Clip Ø mm	Covering
SURECLIP™ SMALL					
ROCC-C-26-195-C	8	1950	135°	2.6	Pink
ROCC-C-26-230-C	8	2300	135°	2.6	Pink
SURECLIP™					
ROCC-D-26-195-C	11	1950	135°	2.6	Blue
ROCC-D-26-230-C	11	2300	135°	2.6	Blue
SURECLIP™ PLUS					
ROCC-F-26-195-C	16	1950	135°	2.6	Yellow
ROCC-F-26-230-C	16	2300	135°	2.6	Yellow
ECO CLIP					
ROCC-D-26-195	11	1950	135°	2.6	Without
ROCC-D-26-230	11	2300	135°	2.6	Without
PREMIUM ECO CLIP BLUE					
ROCC-D-26-195-C-S	11	1950	135°	2.6	Blue
ROCC-D-26-230-C-S	11	2300	135°	2.6	Blue
PREMIUM ECO CLIP YELLOW					
ROCC-F-26-195-C-S	16	1950	135°	2.6	Yellow
ROCC-F-26-230-C-S	16	2300	135°	2.6	Yellow

Packaging unit: 10 pieces

